

User Management

Value	Access Level	Detailed Description
	Complete access	Ability to view content Ability to create/edit content where applicable
	Read-only access	Ability to view content No ability to create/edit content where applicable
	No access	No ability to view content No ability to create/edit content where applicable

User Management

Portal Component

Basic Clinic Access

Dashboard and Reports	Green
Templates	Yellow
Users	Yellow
Patient Snapshot and Treatment Summary	Green
Device Programs, System and Patient Settings	Red
Device Program, System Settings and Patient Settings Confirmation	Yellow
Device Settings List	Red
Device Programs, System and Patient Templates	Red
Device Templates List	Red
Flag rules	Red
Clinic user information, roles, patient access, search and list	Yellow
User creation confirmation	Yellow
Forgot password – enter email, enter security answers, create password and security questions and change password	Green
Account settings	Green
Change username and name	Green
Edit preferences	Green
Login, help, FAQs, Contact Us	Green
Privacy Policy/Terms and conditions	Green

User Management

Portal Component

Homechoice Claria Device Manager

Dashboard and Reports	Yellow
Templates	Yellow
Users	Yellow
Patient Snapshot and Treatment Summary	Yellow
Device Programs, System and Patient Settings	Green
Device Program, System Settings and Patient Settings Confirmation	Green
Device Settings List	Green
Device Programs, System and Patient Templates	Yellow
Device Templates List	Yellow
Flag rules	Yellow
Clinic user information, roles, patient access, search and list	Yellow
User creation confirmation	Yellow
Forgot password – enter email, enter security answers, create password and security questions and change password	Green
Account settings	Green
Change username and name	Green
Edit preferences	Green
Login, help, FAQs, Contact Us	Green
Privacy Policy/Terms and conditions	Green

User Management

Portal Component

Homechoice Claria Clinical Settings Manager

Dashboard and Reports	Yellow
Templates	Green
Users	Yellow
Patient Snapshot and Treatment Summary	Yellow
Device Programs, System and Patient Settings	Yellow
Device Program, System Settings and Patient Settings Confirmation	Yellow
Device Settings List	Yellow
Device Programs, System and Patient Templates	Green
Device Templates List	Green
Flag rules	Green
Clinic user information, roles, patient access, search and list	Yellow
User creation confirmation	Yellow
Forgot password – enter email, enter security answers, create password and security questions and change password	Green
Account settings	Green
Change username and name	Green
Edit preferences	Green
Login, help, FAQs, Contact Us	Green
Privacy Policy/Terms and conditions	Green

User Management

Portal Component

Clinic User Manager

Dashboard and Reports	Yellow
Templates	Yellow
Users	Green
Patient Snapshot and Treatment Summary	Yellow
Device Programs, System and Patient Settings	Yellow
Device Program, System Settings and Patient Settings Confirmation	Yellow
Device Settings List	Yellow
Device Programs, System and Patient Templates	Green
Device Templates List	Green
Flag rules	Green
Clinic user information, roles, patient access, search and list	Yellow
User creation confirmation	Yellow
Forgot password – enter email, enter security answers, create password and security questions and change password	Green
Account settings	Green
Change username and name	Green
Edit preferences	Green
Login, help, FAQs, Contact Us	Green
Privacy Policy/Terms and conditions	Green